

CASSETTE CAPACITY CHART



CLEANING EQUIPMENT CAPACITY LIST

MANUFACTURER/PHONE	MODEL	TANK DIMENSIONS L x W x D (INCHES)	CAPACITY IN GALLONS	IMS® CASSETTES													RESIN CASSETTES	
				EXAM & HANDPIECE 3 x 8 x 1.25	SMALL 10 INSTRUMENT 5.5 x 8 x 1.25	TWO-TIER 8-INSTRUMENT 3 x 8 x 1.5	SPACE SAVER SEVEN 4.5 x 8 x 1	TWO-TIER 14-INSTRUMENT 4.5 x 8 x 1.5	SIGNA STAT 12 INSTRUMENT 6.5 x 10.5 x 1.25	TWO-TIER 20-INSTRUMENT 6 x 9.7 x 1.5	LARGE 8 x 11 x 1.25	ORAL SURGERY 8 x 14.5 x 1.25	ORTHO BOX 7.25 x 5.75 x 3.25	SMALL THIN ORTHO 7 x 8 x 0.75	LARGE THIN ORTHO 8 x 11 x 0.75	TWO-TIER ORTHO 20-INSTRUMENT 6 x 8 x 1.5	SMALL & TWO-TIER 4.5 x 8 x 1.5	LARGE 8 x 11 x 1.5
Brand Max ® 5	U-10LH / U-10LHREC	12 x 9.5 x 6	2.64	12	7	9	10	6	4	3	N/C	N/C	4	6	N/C	3	6	N/C
	U-13LH / U-13LHREC	13 x 11.8 x 6	3.43	12	8	9	10	6	4	3	4	N/C	4	6	6	3	6	3
	U-20LH / U-20LHREC	19.5 x 11.8 x 6	5.28	24	12	18	16	12	8	6	4	4	12	12	6	9	12	3
Biotrol 1-800-822-8550	Purit 300	12.25 x 11.25 x 5.5	3	12	6	12	9	7	3	3	3	N/C	1	5	5	3	7	3
	Purit 150	11 x 5.25 x 5.625	1.5	4	1	4	4	3	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	3	N/C
Coltene Whaledent 1-800-221-3046	UC300 (cassette rack accessory P/N UC320 or UC325)	14 x 9 x 6	3	9	6 [^]	9	6	6	3	3	3	N/C	3	3	3	3	6	3
	UC300R (cassette rack accessory P/N UC320 or UC325)	14 x 9 x 6	3	9	6 [^]	9	6	6	3	3	3	N/C	3	3	3	3	6	3
	UC125 (UC164XD)	11.25 x 5.62 x 5	1.25	1	N/C	1	3	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	1	1	N/C
	UC125 (UC151XD - with original basket) •	11.25 x 5.62 x 5	1.25	2	N/C	1	2	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	1	N/C
	UC125 (UC151XD - with updated basket) ›	11.25 x 5.62 x 5	1.25	4	N/C	4	3	2	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	2	N/C
	UC100	11.5 x 5.25 x 4	1.25	1	N/C	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C
Esma, Inc. 1-800-276-2466	E789	13 x 11 x 11	6.5	20	12	12	12	9	8	8	6	N/C	4	10	10	8	9	6
	E789-oms	15 x 10 x 10	6.5							8	6	6	4	10	10	10	9	6
	E397	12 x 10 x 8	3.5	12	7	9	8	6	5	4	3	N/C	2	7	7	5	6	3
	E289	12 x 10 x 8	3.5	12	7	9	8	6	5	4	3	N/C	2	7	7	4	6	3
	E888R	12 x 10 x 8	3.5	12	7	9	8	6	5	4	3	N/C	2	7	7	4	6	3
	E888	12 x 10 x 8	3.5	12	7	9	8	6	5	4	3	N/C	2	7	7	4	6	3
Health Sonics 1-800-342-3096	T19.9	20.8 x 13.5 x 12.2	5	17	9	17	16	8	3	2	3	3	2	10	5	4	8	3
	T13.7	14 x 14 x 10.7	3.5	12	6	12	9	8	3	2	3	N/C	1	5	5	2	8	3
	T10.4	13 x 11.5 x 10.3	2.7	10	3	10	8	6	N/C	3	N/C	N/C	1	5	N/C	3	6	
	T10.4 C/S	13 x 11.5 x 12.3	3.2	12	7	12	12	8	N/C	4	N/C	N/C	2	8	N/C	4	8	
	T9	20.8 x 7.5 x 11.5	2.3	6	N/C	6	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C
	T5.7	12.8 x 8 x 10.5	1.5	4	N/C	4	4	3	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	3	
	T4.4	10.5 x 7.5 x 10	1.1	1	N/C	1	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C
	T3.3	10.5 x 7.5 x 8.5	0.8	2	N/C	2	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C


All capacity quantities and dimensions are provided by listed equipment manufacturers. Hu-Friedy provides this for your convenience, but we encourage you to verify with equipment manufacturer before purchase. We are happy to provide sample cassettes to confirm fit.

N/C = Not Compatible
 ■ Must be Hu-Friedy part code
 Infinity Series - IMN416 (color) 0 - 8 x 10.5
 Signature Series - IM416 (color) 0 - 8 x 10.5
 to accommodate rounded corner of tank

[^] According to equipment Manufacturer, when the UC300/UC300R is used with Hu-Friedy 10-instrument cassettes, in order to obtain full cassette coverage, it may be necessary to fill the tank beyond the fill line. Ensure the water level does not exceed the lip of the plastic bezel inside the tank. This action will not void the standard warranty.

+ Recessed models
[†] Must order 2 x 4XL Cassette Racks (Full Size) - part #01-109965S
 • Original style basket with integrated instrument holders on bottom
 › New style basket with flat bottom

[§] In order to obtain full cassette coverage, it may be necessary to fill tank beyond the fill line. Please ensure water does not extend above top of tank. Capacity show with optimal instrument cassette basket.



CASSETTE CAPACITY CHART



CLEANING EQUIPMENT CAPACITY LIST (CONTINUED)

MANUFACTURER/PHONE	MODEL	TANK DIMENSIONS L x W x D (INCHES)	CAPACITY IN GALLONS	IMS® CASSETTES													RESIN CASSETTES	
				EXAM & HANDPIECE 3 x 8 x 1.25	SMALL 10 INSTRUMENT 5.5 x 8 x 1.25	TWO-TIER 8-INSTRUMENT 3 x 8 x 1.5	SPACE SAVER SEVEN 4.5 x 8 x 1	TWO-TIER 14-INSTRUMENT 4.5 x 8 x 1.5	SIGNA STAT 12 INSTRUMENT 6.5 x 10.5 x 1.25	TWO-TIER 20-INSTRUMENT 6 x 9.7 x 1.5	LARGE 8 x 11 x 1.25	ORAL SURGERY 8 x 14.5 x 1.25	ORTHO BOX 7.25 x 5.75 x 3.25	SMALL THIN ORTHO 7 x 8 x 0.75	LARGE THIN ORTHO 8 x 11 x 0.75	TWO-TIER ORTHO 20-INSTRUMENT 6 x 8 x 1.5	SMALL & TWO-TIER 4.5 x 8 x 1.5	LARGE 8 x 11 x 1.5
L&R 201-991-5330	SweepZone 200 (No Rack - Use Basket)	11.75 x 6 x 4	1.06	2	N/C	1	1	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	2	N/C
	SweepZone 310R / SweepZone 310 / Quantrex 310 (Rack #26665)*	15.5 x 9 x 6	3.25	12	6	8	6	4	3	2	3	N/C	2	8	4	4	4	2
	SweepZone 360R / Quantrex 360 Rack #19368*	11.5 x 9.25 x 8	3.6	12	4	12	10	8	4	4	N/C	N/C	2	7	N/C	3	6	N/C
	SweepZone 650R / Quantrex 650 Rack #26847 (Holds 2 Racks)*	19.75 x 11.75 x 6.5	6.5	18	6	18	16	12	6	6	6	N/C	2	10	10	6	12	6
Midmark Soniclean 1-800-MIDMARK 1-800-643-6275	Soniclean M550 (rack# 9A407001)	19.5 x 11.5 x 6	5.5	15	9	15	16	12	6	6	6	3	3	10	10	9	12	6
	M550 with basket			15	9	15	16	12	6	6	3	3	3	10	5	9	12	3
	Soniclean M250 (rack# 9A295001)	11.5 x 9.5 x 6	2.5	9	3	9	8	4	3	3	3 [■]	N/C	1	6	N/C	3	6	3 [■]
	M250 with basket			9	3	9	8	4	3	3	0	0	1	6	5	3	6	1
	Soniclean M150	11.5 x 6 x 6	1.5	3	0	3	4	3	N/C	N/C	N/C	N/C	1	N/C	N/C	N/C	3	N/C
Miele Canada 1-888-325-3957 USA 1-80-991-9380	G7881/G7891 (Rack: 4xE523)	19.6 x 21 x 19.6	N/A	55	35	55	35	20	11	55	11	11	12	24	11	24	35	11
Quala 1-800-929-4ADC	Sweep 5300 (Rack# 26665 for large cassettes)	15.5 x 9 x 6	3.25	12	6	8	8	6	3	2	3	N/C	1	5	5	6	6	3
	Sweep 5200	11.75 x 6 x 4	1.06	2	N/C	2	2	0	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	2	N/C
	5002	9.4 x 5.4 x 4	0.84	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C
SciCan 1-800-667-7733	Hydrim L110W G4 (**w/ modified racks)	23 x 15.5 x 18	N/A	30	20	30	15	15	5/**10	5/**10	**10	**8	8	10	**10	10	15	**10
	New Hydrim C61W G4 (Std rack / 6 rack)	12 x 19.5 x 16	N/A	12/24	6/12	12/24	9/18	9/18	3/6	3/6	3/6	3/6	2	6/12	3/6	6/12	9/18	3/6
Sonix IV 1-843-554-0239 1-800-878-2769	SR118/ST118	19.5 x 11.5 x 8	7.8	24	12	20	20	16	4	10	4	4	5	12	12	10	16	4
	SR116/ST116	19.5 x 11.5 x 6	5.5	18	9	15	14	12	3	6	3	3	5	8	8	6	12	3
	SR236/ST236	11.5 x 12.8 x 6	3.5	9	6	9	10	7	3	3	3	N/C	2	4	4	3	7	3
	SR128/ST128	11.5 x 9.5 x 8	3.5	12	8	12	10	8	4	4	4	N/C	2	6	6	6	8	3
	SR126/ST126	11.5 x 9.5 x 6	2.7	6	3	9	6	6	3	3	N/C	N/C	1	4	N/C	3	6	N/C
Sultan 1-800-637-8582	ProSonic 2000 (2 Racks #60230)	19.5 x 11.5 x 6	5.5	18	12	18	16	12	6	6	6	4	4	10	10	6	12	6
	ProSonic 1000 (1 Rack #60230)	11.5 x 9.5 x 6	2.5	12	4	9	8	6	3	3	3 [*]	N/C	2	5	N/C	3	6	3
	ProSonic 600	11.5 x 6 x 6	1.5	4	4	3	4	3	N/C	3	N/C	N/C	1	N/C	N/C	N/C	3	N/C
	ProSonic 300	9.5 x 5.5 x 4	0.75	3	N/C	2	N/C	2	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	2	N/C
Tuttner 1-800-624-5836	CSU3	15.5 x 9 x 6	3	9	6	6	6	6	4	3	3	N/C	4	5	5	6	6	3
	CSU1	9.38 x 5.38 x 4	0.8	2	N/C	2	2	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	1	N/C


All capacity quantities and dimensions are provided by listed equipment manufacturers. Hu-Friedy provides this for your convenience, but we encourage you to verify with equipment manufacturer before purchase. We are happy to provide sample cassettes to confirm fit.

N/C = Not Compatible
 ■ Must be Hu-Friedy part code
 Infinity Series - IMN416 (color) 0 - 8 x 10.5
 Signature Series - IM416 (color) 0 - 8 x 10.5
 to accommodate rounded corner of tank

* According to equipment Manufacturer, when the UC300/UC300R is used with Hu-Friedy 10-instrument cassettes, in order to obtain full cassette coverage, it may be necessary to fill the tank beyond the fill line. Ensure the water level does not exceed the lip of the plastic bezel inside the tank. This action will not void the standard warranty.

+ Recessed models
 † Must order 2 x 4XL Cassette Racks (Full Size) - part #01-109965S
 • Original style basket with integrated instrument holders on bottom
 > New style basket with flat bottom

§ In order to obtain full cassette coverage, it may be necessary to fill tank beyond the fill line. Please ensure water does not extend above top of tank. Capacity show with optimal instrument cassette basket.



CASSETTE CAPACITY CHART



STERILIZER EQUIPMENT CAPACITY LIST

MANUFACTURER/PHONE	MODEL	CHAMBER DIMENSIONS L x W x D (INCHES)	IMS® CASSETTES													RESIN CASSETTES	
			EXAM & HANDPIECE 3 x 8 x 1.25	SMALL 10 INSTRUMENT 5.5 x 8 x 1.25	TWO-TIER 8-INSTRUMENT 3 x 8 x 1.5	SPACE SAVER SEVEN 4.5 x 8 x 1	TWO-TIER 14-INSTRUMENT 4.5 x 8 x 1.5	SIGNA STAT 12 INSTRUMENT 6.5 x 10.5 x 1.25	TWO-TIER 20-INSTRUMENT 6 x 9.7 x 1.5	LARGE 8 x 11 x 1.25	ORAL SURGERY 8 x 14.5 x 1.25	ORTHO BOX 7.25 x 5.75 x 3.25	SMALL THIN ORTHO 7 x 8 x 0.75	LARGE THIN ORTHO 8 x 11 x 0.75	TWO-TIER ORTHO 20-INSTRUMENT 6 x 8 x 1.5	SMALL & TWO-TIER 4.5 x 8 x 1.5	LARGE 8 x 11 x 1.5
A-dec / W&H 1-800-547-1883	LISA MB 17	9.8 x 13.8	9	6	N/C	6	N/C	N/C	N/C	3	N/C	N/C	3	3	N/C	N/C	N/C
	LISA MB 22 (Available in Canada)	9.8 x 17.7	15	6	N/C	9	N/C	3	N/C	3	3	N/C	6	3	N/C	N/C	N/C
Midmark 1-800-MIDMARK 1-800-643-6275	M11 w/ Horizontal Rack #9A215001	11 x 18	16	8	16	12	12	1	1	4	4	2^	8	4	0	12	4
	M11 w/ Vertical Rack #9A215002	11 x 18	16	8	16	12	12	4	4	4	4	N/C	8	4	4	12	4
	M11 with trays	11 x 18	12	6	12	8	8	3	4	2	2	3	6	3	4	8	2
	M9 standard rack with trays	9 x 15	6	3	6	4	3	2	2	1	N/C	2	3	2	2	4	1
	M7 standard rack with tray	7.5 x 15	3	2	3	2	2	N/C	1	N/C	N/C	N/C	N/C	N/C	1	2	N/C
	M3 with tray	7.6 x 1.6 x 12.1	2	1	N/C	1	N/C	1	N/C	N/C	N/C	N/C	1	N/C	N/C	N/C	N/C
Miele Canada 1-888-325-3957 USA 1-80-991-9380	PS1201B, Class B	10x16	9	6	6	6	3	3	3	3	2	N/C	6	3	3	3	N/C
Pelton & Crane 1-800-659-6560	Magnaclave	15.25 x 15 x 26	28	28	28	28	28	14	14	14	3	6	14	14	14	28	14
	Delta XL10 w/ rack (024338)	9.8 x 17.5	15	9	12	9	9	3	3	3	3	2	4	3	9	9	3
	Delta XL8	8.25 x 14	8	1	8	1	N/A	0	0	N/C	N/C	1	2	N/C	1	4	N/C
	Delta Q 10 Rack# (024338)	9.8 x 17.5	15	9	12	9	9	3	3	3	3	2	4	3	9	9	3
	Delta Q 8	8.25 x 14	8	1	8	1	N/A	0	0	N/C	N/C	1	2	N/C	1	4	N/C
SciCan Canada 1-800-667-7733 USA 1-800-572-1211	Bravo17	10 x 13.5	9	6	9	6	6	3□	3	3	N/A	3	3	3	3	6	3
	Bravo17V	10 x 13.5	9	6	9	6	6	3□	3	3	N/A	3	3	3	3	6	3
	Bravo21V	10 x 17.5	12	6	12	9	9	3□	3	3	3	3	6	3	6	9	3
	Quantim 16	9.8 x 13	9	6	9	6	6	3	3	3	N/C	1	3	3	3	6	3
	Statim 7000	8.6 x 13.5 x 2.5	4	4	4	4	2	2	1	2	N/C	N/C	3	3	3	2	2
	Statim 5000 G4	7.1 x 15 x 3.2	2	2	2	2	1	2	1	N/C	N/C	N/C	3	N/C	N/A	1	N/C
	Statim 2000 G4	7.1 x 11 x 1.6	2	1	1	1	1	1	1	N/C	N/C	N/C	1	N/C	N/A	N/C	N/C
Tuttnauer 1-800-624-5836	3870*	15x30	45	30	45	30	30	15**	15**	15	4	10	15**	15**	15**	30	15
	EZ10	10 x 19	16	9	12	9	9	4	4	3	3	2	6	3	6	9	3
	2540M	10 x 19	16	9	12	9	9	4	4	3	3	2	6	3	6	9	3
	EZ9	9 x 18	12	6	8	6	6	3	3	2	2	2	4	2	4	6	2
	2340M	9 x 18	12	6	8	6	6	3	3	2	2	2	4	2	4	6	2
	EZ9Plus	9x19.8	12	6	8	6	6	2	2	2	2	2	4	2	4	6	2
	Valueklave 1730	7 x 12	2	2	2	2	2	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	2	N/C
	EZ11Plus	11x19.8	16	12	16	12	12	4	4	4	4	2	8	4	8	12	4
	Elara 11 Class B	11 x 19.8	20	12	16	12	12	5	5	4	4	2	8	4	8	12	4

All capacity quantities and dimensions are provided by listed equipment manufacturers. Hu-Friedy provides this for your convenience, but we encourage you to verify with equipment manufacturer before purchase. We are happy to provide sample cassettes to confirm fit.

N/C = Not Compatible
 ■ Must be Hu-Friedy part code
 Infinity Series - IMN416 (color) 0 - 8 x 10.5
 Signature Series - IM416 (color) 0 - 8 x 10.5
 to accommodate rounded corner of tank

** Optional Rack required for these Cassettes
 □ Loose on the rails
 ◆ Does not fit within rail
 ^ Insert tray on bottom shelf and then place ortho boxes